

PSA1/PL1 Pillow Loudspeaker System



DESCRIPTION

Designed for application into hospital rooms, the PSA1 selector & attenuator and the PL1 pillow loudspeaker allow that the patients listen to personal audio programs without disturbing their neighbors.

FEATURES

- The PSA1 control plate was designed so as to fit standard hospital equipment bays and uses an inductive attenuator with no easily wearable parts
- 6-program selector
- 6-step attenuator with off position
- ABS enclosure and Mylar speaker on PL1
- PL1 is conveniently small-sized and easily goes through the regular hospital accessories cleaning procedures.
- Fully EN 55103 compliant according to 89/336/EMC

SPECIFICATIONS

Frequency Response	200Hz-8kHz
Nominal Power	0.3W
Connection	1,8m cable, 1/4"(6,35mm) jack
Input	100V screw-terminal type
Attenuation steps	Off, 24,-18, 12, 6, 0dB (Máx.)
PSA1 maximum output	1W / 8 Ohm
Physical	115x115x27mm (PL1) 80x80x50mm (PSA1)
Weight	165g (PL1), 370g (PSA1)
Operating temperature	0-40°C

PSA1/PL1 V1.11 2003-11-10

Specifications may change without prior notice

Copyright © by Acutron 2003 all rights reserved.

LARSEN® is a trademark of Acutron Electroacústica

ABC-Abrunheira Business Center
Edifícios 10-11 R. das Maçarocas
Abrunheira 2710-056 Sintra
Tel +351 219152040
Fax +351219154079
web: www.acutron.net
mail: info@acutron.net

DISTRIBUTOR



LARSEN®